# Liquordiagnostik FSME / TBE Virus

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## CSF analysis FSME / TBE Virus

*recomWell* FSME / TBE Virus IgG and IgM for quantitative determination of intrathecally produced human IgG and IgM antibodies to Tick borne encephalitis (TBE) Virus.

### 1. General aspects / Intended use

In suspected cases of TBE Virus infection, the detection of intrathecally produced FSME / TBE Virus-specific IgG and IgM antibodies is indicated. For this purpose serum / cerebrospinal fluid (CSF) pairs are analysed in parallel in a quantifiable test (e.g. *recomWell* FSME / TBE Virus). An excel program developed to facilitate the calculations for CSF diagnosis (acc. to Reiber) is available from MIKROGEN. See also Chapter 3. For manual calculation a patient protocol is available from MIKROGEN.

First, serum and CSF are included in the standard TBE Virus screening test with *recomWell* FSME / TBE Virus. In case of a positive CSF, the values for total IgG and/or IgM and albumin must be requested for serum and CSF as basis for the determination of the TBE Virus specific antibody index.

By this method based on the FSME / TBE Virus-specific antibody index, an intrathecally produced antibody fraction can be differentiated from a fraction that has passed the blood-CSF barrier passively to enter the subarachnoid space. Applying the limiting value QLim from the quotient diagram according to Reiber (1,2) ensures that no false-negative antibody index values will be calculated in case of a polyspecific intrathecal immune reaction.

### 2. Screening of serum and CSF in *recomWell* FSME / TBE Virus

Serum and CSF must be withdrawn from the patient at the same day. To avoid test-to-test variations, serum and CSF must always be measured in parallel in the same test run. To increase diagnostic safety, patient sera and CSF should be tested in duplicate.

#### 2.1. Sample dilution

<b>IgG</b>	Serum 1:101	In accordance to the instructions for use, the control sera from the kit (positive control, cutoff control, negative control) and the patient sera are diluted 1:101 with dilution buffer. 100 µl of each of these dilutions are pipetted per well.
<b>IgG</b>	CSF 1:4	25 µl CSF + 75 µl dilution buffer  75 µl of dilution buffer are first pipetted into the well, followed by 25 µl of CSF and the solution is then mixed well.

<b>IgM</b>	Serum 1:101	In accordance with the instructions for use, the control sera from the kit (positive control, cutoff control, negative control) and the patient sera are diluted 1:101 with dilution buffer. 100 µl of each of these dilutions are pipetted per well.
<b>IgM</b>	CSF 1:2	50 µl CSF + 50µl dilution buffer  50 µl of dilution buffer are first pipetted into the well, followed by 50 µl of CSF and the solution is then mixed well.

For serum-CSF pairs no RF-absorption is necessary. The further test procedure is performed according to the instructions for use, *recomWell* FSME / TBE Virus IgG/IgM.

## 2.2. Evaluation

The test is valid if the criteria for the control sera are fulfilled according to the instructions for use, *recomWell* FSME / TBE Virus IgG/IgM.

With respect to the qualitative evaluation the marginal values for sera and CSF in IgG and IgM are calculated as following:

<b>IgG / IgM</b>	marginal value for serum:      OD cutoff control
<b>IgG / IgM</b>	marginal value for CSF:            0,5 x OD cutoff control

## 2.3. Test interpretation

### IgG determination:

<b>Result of <i>recomWell</i> FSME / TBE Virus IgG</b>	<b>Evaluation and further action</b>
OD <sub>serum</sub> 1:101 < marginal value for serum and OD <sub>CSF</sub> 1:4 < marginal value for CSF	No detection of TBE Virus spec. IgG antibodies in serum No detection of TBE Virus spec. IgG antibodies in CSF No indication of intrathecally produced antibodies <b>No determination of antibody index AI (IgG)</b>
OD <sub>serum</sub> 1:101 ≥ marginal value for serum and OD <sub>CSF</sub> 1:4 < marginal value for CSF	Detection of TBE Virus spec. IgG antibodies in serum No detection of TBE Virus spec. IgG antibodies in CSF No indication of intrathecally produced antibodies <b>No determination of antibody index AI (IgG)</b>
OD <sub>serum</sub> 1:101 < marginal value for serum and OD <sub>CSF</sub> 1:4 ≥ marginal value for CSF	No detection of TBE Virus spec. IgG antibodies in serum Detection of TBE Virus spec. IgG antibodies in CSF Indication of intrathecally produced antibodies <b>Determination of antibody index AI (IgG)</b>
OD <sub>serum</sub> 1:101 ≥ marginal value for serum and OD <sub>CSF</sub> 1:4 ≥ marginal value for CSF	Detection of TBE Virus spec. IgG antibodies in serum Detection of TBE Virus spec. IgG antibodies in CSF Indication of intrathecally produced antibodies <b>Determination of antibody index AI (IgG)</b>

**IgM determination:**

<b>Result of recomWell FSME / TBE Virus IgM</b>	<b>Evaluation and further action</b>
<p>OD<sub>serum</sub> 1:101 &lt; marginal value for serum</p> <p>and</p> <p>OD<sub>CSF</sub> 1:2 &lt; marginal value for CSF</p>	<p>No detection of TBE Virus spec. IgM antibodies in serum</p> <p>No detection of TBE Virus spec. IgM antibodies in CSF</p> <p>No indication of intrathecally produced antibodies</p> <p><b>No determination of antibody index AI (IgM)</b></p>
<p>OD<sub>serum</sub> 1:101 ≥ marginal value for serum</p> <p>and</p> <p>OD<sub>CSF</sub> 1:2 &lt; marginal value for CSF</p>	<p>Detection of TBE Virus spec. IgM antibodies in serum</p> <p>No detection of TBE Virus spec. IgM antibodies in CSF</p> <p>No indication of intrathecally produced antibodies</p> <p><b>No determination of antibody index AI (IgM)</b></p>
<p>OD<sub>serum</sub> 1:101 &lt; marginal value for serum</p> <p>and</p> <p>OD<sub>CSF</sub> 1:2 ≥ marginal value for CSF</p>	<p>No detection of TBE Virus spec. IgM antibodies in serum</p> <p>Detection of TBE Virus spec. IgM antibodies in CSF</p> <p>Indication of intrathecally produced antibodies</p> <p><b>Determination of antibody index AI (IgM)</b></p>
<p>OD<sub>serum</sub> 1:101 ≥ marginal value for serum</p> <p>and</p> <p>OD<sub>CSF</sub> 1:2 ≥ marginal value for CSF</p>	<p>Detection of TBE Virus spec. IgM antibodies in serum</p> <p>Detection of TBE Virus spec. IgM antibodies in CSF</p> <p>Indication of intrathecally produced antibodies</p> <p><b>Determination of antibody index AI (IgM)</b></p>

**The test can be terminated if serum and CSF give negative IgG and IgM results in screening. The FSME / TBE Virus-specific antibody index cannot be determined.**

**In case the CSF OD exceeds the marginal value, the procedure continues. In the following, IgG and IgM tests will be described in parallel. However, the antibody index is to be determined only for the immunoglobulin class in which a positive screening result for the CSF has been obtained.**

## 2.4. Determination of TBE Virus-specific antibody concentrations in serum and CSF

For the conversion of the ODs of CSF and serum into units, extinction values between 0.100 and 1.8 are required (linear part of a dilution curve).

If the CSF is positive and the value for the serum < 0.100, an OD of 0.100 is assumed in the unit calculation for the serum.

If the OD value for the serum or CSF is higher than 1.8 or even outside the measurable range, the serum-CSF pair must be tested again in a parallel run with an appropriate higher dilution. The following dilutions are recommended:

*recomWell* FSME / TBE Virus IgG

OD	Serum dilution	CSF dilution	
		Preparation in the well*	
Standard dilution	1:101	1:4	75 µl DB + 25 µl CSF
Results of standard dilution (OD):	Recommendation for further dilutions:		
<b>OD</b> ≤ 2,000	No further dilution		
2,000 < <b>OD</b> ≤ 3,000	1:500	1:10 1:20	90 µl DB + 10 µl CSF 95 µl DB + 5 µl CSF
<b>OD</b> > 3,000	1:1000, 1:2000, 1:4000 geometric dilution	Preparation in a tube*	
		1:40, 1:200, 1:1000 geometric dilution	

*recomWell* FSME / TBE Virus IgM

OD	Serum dilution	CSF dilution	
		Preparation in the well*	
Standard dilution	1:101	1:2	50 µl DB + 50 µl CSF
Results of standard dilution (OD):	Recommendation for further dilutions:		
<b>OD</b> ≤ 2,000	No further dilution		
2,000 < <b>OD</b> ≤ 3,000	1:500	1:5	80 µl DB + 20 µl CSF
<b>OD</b> > 3,000	1:1000, 1:2000, 1:4000 geometric dilution	Preparation in a tube*	
		1:10, 1:20, 1:40 geometric dilution	

\* Recommendation: The low dilutions of CSF can be carried out directly in the well on the microtiter plate. The high dilutions of CSF should be carried out in a tube (Polypropylen = PP).

DB: dilution buffer

In practice, the method of one-point quantification using the cutoff control has been proven in *recomWell* FSME / TBE Virus. Using one-point quantification, the ODs of serum and CSF obtained in the evaluable range between 0.100 and 1.8 are converted into units using formula [1]. The formula is valid for IgG and IgM. If, after geometric dilution of serum or CSF, several values are in this range, the OD value closest to 1.000 is selected (highest precision with respect to dilution errors and linearity of the standard curve).

**In this formula the value of the cutoff control (cutoff for serum) of the corresponding IgG or IgM determination must be used as cutoff for serum and for CSF. The marginal value for CSF from the qualitative evaluation is not valid here.**

All calculations following can also be performed by the use of the Excel program; (see chapter 3)

$$[1] \quad \text{U/ml serum or CSF} = \frac{\text{OD}_{\text{Serum resp. CSF}} \times 20 \times \text{dilution factor}^{1)}}{\text{Cutoff}^{2)} \times 101^{3)}}$$

<sup>1)</sup> Dilution factor = reciprocal value of dilution, e.g., dilution factor 500 corresponds to a dilution of 1:500

<sup>2)</sup> For the cutoff the value of the respective cutoff control has to be used for **serum and for CSF**.

<sup>3)</sup> Standard dilution in *recomWell* FSME / TBE Virus IgG/IgM

This formula gives the TBE Virus specific IgG and IgM concentration in serum and CSF in units/ml.

Measurements compared with a standard curve show only slight differences which are not meaningful. If quantification via a standard curve is still required, a geometric dilution series of the positive control is recommended. This requires at least 4 values: 1:100, 1:200, 1:400 and 1:800. The units 100, 50, 25 and 12.5 are arbitrarily assigned to these dilutions and the units of the patient samples are read off on the corresponding curve. The respective dilutions of serum and CSF have to be considered.

Thus, the TBE Virus specific IgG quotient  $Q_{\text{spec}}(\text{IgG})$  and the TBE Virus specific IgM quotient  $Q_{\text{spec}}(\text{IgM})$  are calculated according to Formulas [2] and [3]:

$$[2] \quad \mathbf{Q_{spec} (IgG)} = \text{TBE Virus spec. IgG CSF (U/ml)} / \text{TBE Virus spec. IgG serum (U/ml)}$$

$$[3] \quad \mathbf{Q_{spec} (IgM)} = \text{TBE Virus spec. IgM CSF (U/ml)} / \text{TBE Virus spec. IgM serum (U/ml)}$$

## 2.5. Calculation of $Q_{\text{IgG}}$ , $Q_{\text{IgM}}$ , $Q_{\text{Alb}}$ , $Q_{\text{Lim}}(\text{IgG})$ and $Q_{\text{Lim}}(\text{IgM})$

Total IgG, total IgM and albumin in both serum and CSF are determined by using a suitable method.

From the values obtained, the IgG quotient  $Q_{\text{IgG}}$ , the IgM quotient  $Q_{\text{IgM}}$  and the albumin quotient  $Q_{\text{Alb}}$  are calculated according to Formulas [4], [5], and [6].

$$[4] \quad \mathbf{Q_{IgG}} = \text{IgG-conc. CSF} / \text{IgG-conc. serum}$$

$$[5] \quad \mathbf{Q_{IgM}} = \text{IgM-conc. CSF} / \text{IgM-conc. serum}$$

$$[6] \quad \mathbf{Q_{Alb}} = \text{albumin-conc. CSF} / \text{albumin-conc. serum}$$

From  $Q_{\text{Alb}}$ , the limiting values  $Q_{\text{Lim}}(\text{IgG})$  and  $Q_{\text{Lim}}(\text{IgM})$  are calculated according to Formulas [7] and [8].

$$[7] \quad \mathbf{Q_{Lim} (IgG)} = 0,93 \times \sqrt{Q_{\text{Alb}}^2 + (6 \times 10^{-6})} - 1,7 \times 10^{-3}$$

$$[8] \quad \mathbf{Q_{Lim} (IgM)} = 0,67 \times \sqrt{Q_{\text{Alb}}^2 + (120 \times 10^{-6})} - 7,1 \times 10^{-3}$$

## 2.6. Calculation of the TBE Virus specific IgG antibody index AI (IgG) and the TBE Virus specific IgM antibody index AI (IgM)

The values obtained for QIgG and QLim (IgG) are compared and AI (IgG) is calculated using formula [9] or [10].

Correspondingly, QIgM and QLim (IgM) are compared and AI (IgM) is calculated using formula [11] or [12].

if QIgG < QLim (IgG) [9] <b>AI (IgG) = Qspec (IgG) / QIgG</b>	if QIgM < QLim (IgM) [11] <b>AI (IgM) = Qspec (IgM) / QIgM</b>
if QIgG > QLim (IgG) [10] <b>AI (IgG) = Qspec (IgG) / QLim (IgG)</b>	if QIgM > QLim (IgM) [12] <b>AI (IgM) = Qspec (IgM) / QLim (IgM)</b>

## 2.7. Interpretation of the antibody index values

$0,6 < \mathbf{AI (IgG)} \leq 1,3$	normal	$0,6 < \mathbf{AI (IgM)} \leq 1,3$	normal
$1,3 < \mathbf{AI (IgG)} < 1,5$	borderline	$1,3 < \mathbf{AI (IgM)} < 1,5$	borderline
$\mathbf{AI (IgG)} \geq 1,5$	pathological	$\mathbf{AI (IgM)} \geq 1,5$	pathological

Antibody index values  $0,6 < AI \leq 1,3$  are in the normal range. An intrathecal synthesis of antibodies against TBE Virus is unlikely. The detected TBE Virus antibodies in the CSF originate from the serum and have passively passed into the CSF space through the blood-CSF barrier.

Antibody index values  $<0,6$  indicate an analytical fault (e.g. clinical chemistry). It is recommended to check if serum and CSF were taken on the same day and if sample-mix-up can be excluded. A repetitive determination of the values of clinical chemistry and/or of the TBE Virus specific antibodies might be advisable.

Antibody index values between 1.3 and 1.5 are borderline cases. The serum-CSF pair should be retested. Or if a clinical doubt exists, it is recommended to test another serum-CSF pair in the course of the infection.

Antibody index values  $\geq 1.5$  are to be classified as pathological. They indicate intrathecal synthesis of TBE Virus antibodies.

### 3. Description of the Excel program

Into the first window the patient data, the values of the Clinical Chemistry and the results of the *recomWell* are entered.

Regarding the values from Clinical Chemistry the respective units have to be selected (g/l; mg/l; g/dl; mg/dl)

The results of the *recomWell* may be entered with the corresponding dilution factors. The default value for the serum dilution is 1:101, for the CSF dilution in IgG 1:4, respectively in IgM 1:2. The extinction of the cutoff control is taken over by the program for the CSF calculation as serum and CSF are to be examined in the same run.

From the values of clinical chemistry and the data of the *recomWell* the program automatically calculates the quotients and the antibody indices (AI).

**CSF diagnostics *recomWell*** Version 9.2

***recomWell***

Borrelia

FSME

Lot IgG  IgM

Date

Operator

*recom* Blot Borrelia  
*recom* Line Borrelia  
Calculation of strip dilution

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We will be pleased to send you further literature on serum / CSF diagnosis.

**Liquordiagnostik *recomWell* FSME / TBE Virus IgG, IgM**  
**CSF analysis *recomWell* FSME / TBE Virus IgG, IgM**

Gebrauchsinformation Version/ Instructions for use version:  
gültig ab/ valid from:

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**MIKROGEN** GmbH