


# Fast & easy determination of the immune status for Toxo, Rubella, CMV and HSV-1/2



## On the same strip

- determination of the time of infection (Toxo, CMV)
- protective/non-protective immunity (Rubella)
- HSV-type differentiation (HSV-2)

Toxo lysate	determination of the immune status
Toxo p30	determination of the time of infection

Rubella WHO Cut off control	control for the differentiation of protective and non-protective immunity
Rubella lysate	determination of the immune status

CMV lysate	determination of the immune status
CMV gB2	determination of the time of infection

HSV-1/2 lysate	determination of the immune status
HSV-2 gG2	determination of the HSV-type 2

## Fast & easy determination of IgG and IgM immune status

